




A. Features

- ➤ Wireless receiving heartbeat signal transmitted from heartbeat brooch, chest band frequency is AM5.3KHz;
- ➤ Time display: Hour, Minute, Second;
- Alarm clock function;
- \triangleright Heart rate exceeds the upper limit alarm function ($60 \sim 220$).
- > Timing (maximum time-keeping is 23 hours, 59mimutes, 59 seconds, 99)
- ➤ 12/24 hour displaying format
- ➤ EL backlight
- Keystroke tone

B. Product Function Mode



1. Press key A in any mode, EL background will light for 3 seconds.


2. Time Mode

- In normal time mode, press key B to enter alarm mode; press key D to enter heart monitoring mode; Press key D again to stop timing and display average heart rate and timed time 2 seconds each alternately. In timing stopped mode, press key D for 2 seconds to clear the timed time and heart rate value, then it will return to normal time mode automatically.
- ➤ Time Setting: In time mode, press and hold key C for 3 seconds to access to upper and lower limits setting of time and heart rate (heart rate range is 60~220).

2013-07-04 Model: BJ2954 Page: 1 / 4 Edition: B








- When numerical value flashes, press key D to adjust the flashing value;
- In normal time setting mode, it will exit and return normal time mode automatically if without any operation within 60 seconds.

3. Alarm Mode

- > In alarm mode, press key D to enable/disable alarm () and Hour Announcement () symbol in turn.
- When the set time is reached, the alarm tone will stop automatically after ringing for 30 seconds or stop by pressing any key. "•••) " symbol flashes when alarm tone is heard.
- Alarm Setting: In alarm mode, press key C once to enter alarm setting mode. First, "Hour" flashes. Press key D to adjust the hour; then press key C to flash "minute", and press key D to adjust the minute; finally, press key C to exit setting mode.
- ➤ In alarm setting mode, it will exit and return normal time mode automatically if without any key operation within 60 seconds.

4. Heart Monitoring Mode



- In normal time mode, press key D once to start timing, and at the same time, heart monitoring starts. There will be "BB" alarm tone when out-of-limit heart rate is measured. During timing process, press key D once to stop timing. It displays average heart rate value and timed time for 2 seconds each alternately.
- In heart monitoring mode, when heart rate signal is received, it updates the heart rate value about every 5 seconds; Without heart rate signal, it enters power saving mode 5 seconds later, and heart rate value displays as zero, heart rate symbol disappears.
- In timing stop mode, press key D to continue with timing.
- In timing stop mode, press and hold key D for 2 seconds, the timed time and heart rate value are reset to zero. Then it exits and returns to normal time mode automatically.


5. Use Instruction on Heartbeat Chest Band:

2013-07-04 Model: BJ2954 Page: 2 / 4 Edition : B





- To wear chest band, face the battery cover of the chest band emitter down, then adjust the emitter at center of your chest.
- The chest band is used to record heart beat number and transfer the data to the wrist watch. Please pay attention to the following when wear your chest band:
 - 1. Get the conductive mat on inner side of the chest band wet with several drops of water or conductive gel to ensure good contact.
 - 2. Wear the chest band alone your chest. To ensure accurate heart rate signal, please adjust the binding belt properly and make you chest comfortable.




► Heart rate signal can be transferred within 2 seconds after chest band being worn. Now you can receive heart rate value in heart monitoring mode. "♥" flashes when heartbeat signal is received.

Tips:

- The location of chest band may affect its performance. Please move the chest band alone the binding belt to make it locate above your heart.
- Avoid wear your chest band at the position with bushy sternal seta.
- In dry and cold weather, it may take several seconds for chest band to perform normally. This is normal phenomenon. After several seconds' heating up, it will perform better.
- ➤ As illustrated below, you can either wear your wrist watch on your wrist or fix it on bicycle or exercise machine with a distance ≤65cm from chest band emitter.





Troubleshooting:


If no heart rate value can be received 2 minutes later after chest band being worn, it may be insufficient battery voltage or the dry and cold environment.

2013-07-04 Model: BJ2954 Page: 3 / 4 Edition: B





Solution: ① try replacing battery.



② Wipe your chest with drinking water or conductive gel.

C. Quality Standard

- ➤ Movement size
- **➤** Thickness
- ➤ Work temperature
- ➤ Work voltage
- Veracity
- > Battery type
- > Static maximum drive current
- > Static average drive current
- ➤ Heartbeat monitoring current
- > EL lighting current
- ➤ Alarm ring current
- Battery Life term (El lighting twice per day, alarm ringing once every day for 30 seconds each time; Heart rate measuring for 1hour.)

- : Ф33.2 mm
- : $9.35 \pm 0.15 \text{ mm}$
- : -10°C ~+60°C
- : 3V
- ± 90 seconds/month
- : CR2032
- : ≤3.5 u A (T=25°C)
- : \leq 2.5 u A (T=25 °C)
- : ≤30 u A (T=25 $^{\circ}$ C)
- : \leq 10.0 m A
- $\leq 2.5 \text{ m A}$
- : 36 months (Japan battery)

Attachment: Chest Band

- Quiescent Current
- Emission Transient Current
- > Emission Average Current
- > Emission Distance
- Battery life term (Calculate by using one hour per day.)
- : \leq 1.0 u A (T=25°C)
- : \leq 4.5 m A (T=25 °C)
- $< \le 0.4 \text{ m A} \text{ (T=25°C)}$
- : ≥65 cm
- : 14 months (Japan Battery)

	Prepared	Checked	Approve
Name	冯碧颜	赵星亮	宋检望
Date	13/07/04	13/07/04	13/07/04

2013-07-04 Model: BJ2954 Page: 4 / 4 Edition: B